

QUALITATIVE INORGANIC ANALYSIS

Preliminary Tests

| # COLOUR OF THE SOLID COMPOUND : | INFERENCE |
|----------------------------------|---|
| Blue or bluish green | Cu^{2+} |
| Green | Ni^{2+} |
| Dark brown | Fe^{3+} |
| Light pink or flesh colour | Mn^{2+} |
| Rose red or violet | Co^{2+} |
| White | Absence of $\text{Cu}^{2+}, \text{Ni}^{2+}, \text{Fe}^{3+}, \text{Ca}^{2+}$ |

| # ODOUR : | INFERENCE |
|----------------------------|---------------------------|
| Vinegar link smell | CH_3COO^- |
| Rotten eggs smell | S^{2-} |
| Smell of SO_2 gas | SO_3^{2-} |
| Ammoniacal smell | NH_4^+ |

| # DRY HEATING TEST : | INFERENCE |
|--|--|
| Colorless & odourless gas | CO_3^{2-} |
| Colorless gas with rotten eggs smell | S^{2-} |
| White fumes with suffocating odour | SO_3^{2-} |
| Brown color gas | $\text{NO}^{2-}, \text{NO}^{3-}, \text{Br}^-$ |
| Yellowish green gas with pungent smell | Cl^- |
| Violet vapour | I^- |
| Characteristic vinegar like smell | CH_3COO^- |
| White when hot and blue when cold | Cu^{2+} |
| Yellow when hot and white when cold | Zn^{2+} |
| Blue when hot and pink in cold | CO^{2+} |
| White residue, glows on heating | $\text{Mg}^{2+}, \text{Al}^{3+}$ |
| Black residue | $\text{CO}^{2+}, \text{Fe}^{3+}, \text{Ni}^{2+}$ or Mn^{2+} |

| # COBALT NITRATE TEST: | Inference |
|------------------------|--|
| Green mass | Zn^{2+} |
| Pink mass | Mg^{2+} |
| Blue mass | $\text{Al}^{3+}, \text{PO}_4^{3-}, \text{BO}_3^{3-}$ |
| Dirty green | Sn^{2+} |
| Black | No definite indication |

| # FLAME TEST : | | |
|---|--|------------------|
| Color of flame observed with naked eyes | Color of flame observed through blue glass | INFERENCE |
| Persistent golden yellow | Invisible pink | Na^+ |
| Violet flame of short duration | Pink | K^+ |
| Brick red flame of short duration | Light yellow | Ca^{2+} |
| Apple green flame | Bluish green | Ba^{2+} |
| Crimson | Purple | Sr^{2+} |
| Intense blue surrounded by green flame | Same color visible | Cu^{2+} |

| # Borax Bead Test : Only for coloured salts (copper, iron, manganese, nickel, cobalt, chromium salts) | | | | |
|---|------------|--------------------------------------|---------------|------------------|
| Color of the bead in oxidizing flame | | Color of the bead in reducing flames | | INFERENCE |
| WHEN HOT | WHEN COLD | WHEN HOT | WHEN COLD | |
| Green | blue | Colorless | Red | Cu^{2+} |
| Yellow | Yellow | Green | Yellow | Fe^{3+} |
| Green | Green | Green | Green | Cr^{3+} |
| Pink | Pink | Colorless | Colorless | Mn^{2+} |
| Brown | Pale brown | Grey or black | Grey or black | Ni^{2+} |

ANALYSIS FOR ANIONS

| # DILUTE SULPHURIC ACID TEST: | INFERENCE |
|---|---------------------------|
| Effervescence with colorless & odourless gas turns lime water milk | CO_3^{2-} |
| Brown colored gas turning starch iodide paper blue | NO_2 |
| Colorless gas with rotten egg smell turning lead acetate paper black or grey | S^{2-} |
| Colorless gas with suffocating odour turning potassium dichromate paper green | SO_3^{2-} |
| Colorless vapour having smell of vinegar | CH_3COO^- |

| # POTASSIUM PERMANGANATE TEST: | INFERENCE |
|---|---------------------------|
| Purple colour of KMnO_4 solution is discharged | |
| Yellowish green gas evolved | Cl^- (may be) |
| Red brown gas evolved | Br^- (may be) |
| Violet vapours | I^- (may be) |
| No gas evolved but purple colour is discharged | Fe^{2+} (may be) |

| # CONCENTRATED SULPHURIC ACID TEST: | INFERENCE |
|---|-----------------|
| Colorless gas pungent small. Gives dense white fumes with rod dipped in NH_4OH | Cl^- |
| Brownish gas with pungent smell. Intensity of fumes increases on adding solid MnO_2 | Br^- |
| Evolution of violet vapours turning starch paper blue | I^- |
| Brown fumes intensifying on adding Cu turning or bits of filter paper | NO_3^- |

CONFIRMATORY TEST OF ANIONS

| TEST | OBSERVATION |
|--|--|
| 1. CARBONATE (CO_3^{2-}) # 0.2g salts + 2ml Dil. H_2SO_4 pass through lime water # 1ml W.E + equal amount of MgSO_4 | Lime water turns milky due to formation of CaCO_3 White precipitates of MgCO_3 |
| 2. SULPHITE (SO_3^{2-}) # About 1ml of W.E. + Acetic acid + 2ml BaCl_2 # 1ml of PPT. solution obtained + 2 drops of KMnO_4 solution. | White ppt. soluble in dil. HCl Pink color of KMnO_4 discharges |
| 3. SULPHIDE (S^{2-}) # 1ml of W.E. + a few drops of freshly prepared sod. Nitroprusside solution. # 1ml W.E. + Acetic acid + 1ml lead acetate solution. # shake about 2ml of S.E. with about 0.2g of freshly prepared cadmium carbonate | Purple or violet coloration Black ppt. of PbS . Formation of yellow ppt. of CdS . |
| 4. NITRITE (NO_2^-) # 2ml of W.E. + dil. Acetic acid + a pinch of thiourea and few drops of ferric chloride sol. # 3-4 drop of aq. Sol. + 2 drop of sulphanilic acid + 2 drop of naphthylamine into it # 1ml of W. Extract + few drop of diphenylamine reagent solution. | A blood red coloration Pink coloration A deep blue coloration |
| 5. ACETATE (CH_3COO^-) # 1 ml W. Extract + 1ml of neutral ferric chloride # Heat the pinch of mixture with few drops of conc. H_2SO_4 and 0.5ml of ethyl alcohol in a test tube. | A blood red coloration Pleasant fruity smell due to the formation of ethyl acetate. |
| 6. CHLORIDE (Cl^-) # 1ml of S.E. + dil. HNO_3 + 1ml of AgNO_3 sol. In it. # Heat about 0.2g of the solid mixture with 0.4g of MnO_2 + 2ml conc. H_2SO_4 . # 0.2g of mixture + 0.3g potassium dichromate + 3ml conc. H_2SO_4 and heat pass reddish gas through NaOH sol. + acetic acid + few drop of lead acetate sol. In it. | Crudy white ppt. soluble in NH_4OH Evolution of greenish yellow gas With pungent smell. Yellow ppt. of lead chromate. |
| 7. BROMIDE (Br^-) # 2ml W.E. + dil. HNO_3 + 2ml of AgNO_3 solution. # 2ml of W.E. + Acetic acid + 1-2ml of CS or CHCl_3 to it. Add freshly prepared Cl_2 water or dil. HNO_3 with constant shaking in the tube. | Light yellow ppt. partially soluble in NH_4OH A brown color in the organic layer. |
| 8. IODIDE (I^-) # Neutralize 1-2 ml of the S.E. with acidic acid + few drops of silver nitrate solution. # 2ml of S.E. + dil. Acetic acid + 2ml of CS or CCl_4 drop wise with shaking Cl_2 water or HNO_3 . | Pale yellow ppt. insoluble in NH_4OH A violet color in the organic layer confirms iodide. |
| 9. NITRATE (NO_3^-) # 1ml of W.E. + diphenylamine reagent solution. # 2-3ml of W.E. + acetic acid + 10 drops of ferrous sulphate + 2ml of conc. H_2SO_4 . | Blue color. Formation of brown ring at the junction of two layers. |
| 10. SULPHATE (SO_4^{2-}) # 1ml of W.E. + dil. HCl + 2ml of BaCl_2 solution in it. # 2ml W.E. + acetic acid + few drops of lead acetate solution in it. | White ppt. insoluble in conc. HNO_3 A white ppt. soluble in excess of hot ammonium acetate solution. |
| 11. PHOSPHATE (PO_4^{3-}) # 2ml W.E. + dil. HNO_3 + 1ml of conc. HNO_3 + pinch of solid ammonium molybdate into it. # 2ml W.E. + dil HCl boil off CO_2 + few drops of magnesia mixture solution. | Yellow ppt. White crystalline ppt. |

ANALYSIS FOR CATIONS

To the solution of compound add dil. HCL

| | | | | |
|---------------------------------------|--|---|--|--|
| White ppt. [Group I] Pb ²⁺ | Filtrate | Pass H ₂ S gas through the filtrate | | |
| | ppt. [Group II] Black: Pb ²⁺ , Cu ²⁺ , Hg ²⁺ Yellow: Cd ²⁺ , As ³⁺ | Filtrate. Boil off H ₂ S. Add 3-4 drops of conc. HNO ₃ +add solid ammonium chloride. Boil again, cool & add ammonium hydroxide until the solution becomes basic. | | |
| | ppt. [Group III] Reddish brown: Fe ³⁺ White: Al ³⁺ | Filtrate. Add solid NH ₄ Cl+NH ₄ OH+pass H ₂ S gas | | |
| | ppt. [Group IV] Black: Ni ²⁺ , CO ²⁺ Flesh colored: Mn ²⁺ Dirty white Zn ²⁺ | Filtrate. Boil off H ₂ S+add NH ₄ Cl & NH ₄ OH+(NH ₄) ₂ CO ₃ solution warm and allow to stand. | | |
| | White PPT. [Group V] Ba ²⁺ , Sr ⁺ & Ca ²⁺ | Filtrate. Add NH ₄ Cl, NH ₄ OH & ammonium hydrogen phosphate. White PPT.: Mg ²⁺ Filtrate. Heat to dryness & test for Na ⁺ & K ⁺ with the residue. Original mixture. Heat with NaOH solution. Evolution NH ₃ , NH ₄ ⁺ | | |

ANALYSIS OF GROUP ZERO (NH₄⁺)

| Experiment | Observation | Inference |
|-------------------|--------------------|--|
| NaOH Test | Dense white fumes | NH ₄ ⁺ confirmed |
| Nessler's reagent | Reddish brown ppt. | NH ₄ ⁺ confirmed |

ANALYSIS OF GROUP I [Pb²⁺]

Boil the white precipitate with 5.10ml of water. Divide the solution into three parts.

Confirmation

1. Cool one part under tap. White crystalline ppt. separate out.
2. Potassium iodide test to the second part, add KI solution-yellow ppt. confirms Pb²⁺.
3. Potassium chromate test to the third part add K₂CrO₄ solution -yellow ppt. confirms Pb²⁺.

ANALYSIS OF BLACK PRECIPITATE (PbS, CuS, HgS)

Group II cations can be put into two groups on the basis of color of their sulphates.

| BLACK RESIDUE | FILTRATE | |
|---|---|--|
| | Bluish green | Colorless |
| HgS ppt. Doesn't dissolve in 50% hot nitric acid. So if black residue is left, test Hg ²⁺ ions. | Test for Cu ²⁺ ions. Divide the solution into two parts | Test for Pb ²⁺ ions. To the sol. Add 1ml of ethyl alcohol+2ml conc. H ₂ SO ₄ . White ppt. indicates presence of Pb ²⁺ . Dissolve the white ppt. in 4-5 ml of hot saturated ammonium acetate sol. & divide the solution |
| 1. To one part adds stannous chloride. A white ppt. Turing gray confirms Hg ²⁺ ion presence. | 1.To one part add NH ₄ OH solution in excess deep blue color confirms Cu ²⁺ . | 1. To one part, add potassium chromate solution a yellow ppt. confirms Pb ²⁺ ions. |
| 2. To the other part, add a piece of clean copper wire or foil and set aside. A gray deposit, which on polishing with filter paper becomes silvery white, confirms Hg presence. | 2. Acidify the second part with acetic acid +potassium ferrocyanide solution chocolate color PPT. confirms Cu ²⁺ presence. | 2.to second part add KI solution. Formation of yellow ppt. confirms Pb ²⁺ ions. |

Identification of Cd²⁺ & As³⁺ [yellow ppt.] add about 3-4 ml of yellow ammonium sulphide to ppt. and warm the test tubes for 3 minutes.

| CADMIUM [Cd ²⁺] | ARSENIC [As ³⁺] |
|--|---|
| Reside indicates the presence of Cd ²⁺ dissolve the ppt. In dill. HCl and divide solution into two parts. | Dissolution of the ppt. Indicates presence of As ³⁺ . Add dil. HCl to above solution. Appearance of yellow ppt. in hot conc. HNO ₃ and divide the solution in two parts. |
| i) Dilute one part about 5 times with distilled water+pass H ₂ S gas through the resultant sol. Yellow ppt. confirms Cd ²⁺ . | i) To one part+ammonium molybdate sol. Formation of yellow ppt. confirms As ³⁺ . |
| ii) To the second part add sodium hydroxide sol. Formation of white ppt. Confirms Cd ²⁺ . | ii) To the second part add ammonium hydroxide + magnesia mixture. Appearance of white ppt. confirms As ³⁺ . |

ANALYSIS OF GROUP III [CATIONS]

| Reddish brown ppt. –Fe ³⁺ present | White gelatin ppt. –Al ³⁺ present |
|---|---|
| Dissolve the ppt. In minimum amount of HCl and dissolve the sol. In two parts | Dissolve the ppt. By boiling with NaOH. divide the solution in two parts |
| i) To one part add potassium ferrocyanide solution. A blue color confirms Fe ³⁺ . | i) To one part add NH ₄ Cl and boil. Formation of white gelatin ppt. confirms Al ³⁺ . |
| ii) To the second part add potassium sulphocyanide solution. Appearance of deep red coloration confirms Fe ³⁺ . | ii) Take second part, acidify with dil. HCl. Add two drops of blue litmus followed by ammonia till it smell ammonia. Blue ppt. Floating in colorless solution confirms Al ³⁺ . |
| <p><i>Distinction between Fe²⁺ and Fe³⁺. Add a few drops of potassium ferricyanide solution to the O.S. before adding the conc. HNO₃.</i></p> <p><i>Appearance of deep blue colour or ppt. indicates the presence of Fe²⁺.</i></p> <p><i>Appearance of a deep blue coloration or ppt. On adding potassium ferrocyanide sol. indicates the presence of Fe³⁺ to the original solution.</i></p> | |

ANALYSIS OF GROUP IV [CATIONS]

| BLACK RESIDUE | Pink / White ppt. Dis. ppt in HCl | |
|--|---|---|
| Dissolve the minimum quantity of aqua regia. Evaporate the solution to almost dryness in a porcelain dish. Extract the residue with 5ml of distilled water. Change in color from blue to pink confirms CO ²⁺ . Divide the extract in two parts. | Boil to expel H ₂ S gas. Cool it and add excess of NaOH solution into it. | |
| i) To one part add 1ml of ethyl alcohol+100gm of NH ₄ SCN. shake well. Blue color in alcohol layer confirms the presence of CO ²⁺ in the mixture. | Dissolve the residue in conc. HNO ₃ + 0.2g PbO and boil. Purple solution confirms Mn ²⁺ . | Pass H ₂ S gas through it. White ppt. Confirms Zn ²⁺ presence in the mixture. |
| ii) To the second part add 5ml of dimethyl glyoxime and NH ₄ OH. Pink ppt. confirms the presence of Ni ²⁺ in the mixture. | | |

ANALYSIS OF GROUP V [CATIONS]

| RESIDUE | FILTRATE | |
|--|---|---|
| Residue BaCrO ₄ , wash the ppt. with hot water and perform the flam test. Apple green flame confirms the presence of Ba ²⁺ . | Filter may contain Sr ²⁺ & Ca ²⁺ . Add excess of ammonium sulphate solution & heat. Allow it to stand & filter. | |
| | Residue white SrSO ₄ confirms by flame test. Persistent crimson red flame indicates Sr ²⁺ . | Filtrate. Add ammonium oxalate solution & warm. White ppt. dissolves in HCl & confirms by flame test. Brick red flame Ca ²⁺ confirmed. |

IDENTIFICATION OF Mg²⁺

Divide the filtrate of Group V in two parts.

To one part add 1-2ml of disodium hydrogen phosphate solution followed by ammonia solution till mixture smells ammonia – White crystalline precipitate .

To the second part add magneson reagent (an alkaline solution of p-nitrobenzene azo resorcinol) and excess of NaOH solution – Sky blue precipitate.

Some Hints for Analysing a Salt :

- *If the salt is white in colour ,then Copper ,iron, cobalt , nickel cations are not present.*
- *If the salt contains sulphate as acid radical (anion) then Ba²⁺, Sr²⁺, Ca²⁺ and Pb²⁺ should be absent. (Sulphates of these cations are insoluble & insoluble salts are not in the syllabus.*
- *When Phosphate is confirmed as anion , radicals of Group IV & V & Magnesium should be absent.*
- *To ensure the presence of Group V , take about 1 ml of Group II filtrate & add about 1 ml alcohol & then slowly add 1ml of H₂SO₄.*

Appearance of white precipitates indicates the presence of Group V.
